

Checklist for return of GalliaPharm® in the US

Fill one checklist per package and send one copy to galliapharm@ezag.com for prior-check and before shipping the generator!

Support for completion is given in the document "Explanation of return procedure for GalliaPharm®"

Package (please check and confirm):

- Confirmation that column was completely emptied from liquid
- Confirmation that inlet & outlet have been closed with respective plugs
- Confirmation that GalliaPharm® is placed in proper & qualified plastic bag ("pressure bag"), before packing into the Type A package
- Eckert & Ziegler original Type A package is used
- Adhere all labels on the Type A package (see also the Sheet 3 & Return Authorization Label of your Return Kit)
- A copy of the Type-A container Certificate is part of the package
- A copy of the test report for the Type-A container is part of the package

Check before shipment

Visual inspection of each single part: package, inserts, bag etc. ok

Visual inspection of assembled Type A package ok

Shipment

Shipper/Customer:		Enter EZR RA No.:	
Date:		Current activity at date:	MBq
Serial / Lot No. GalliaPharm®:		Quantity: 1	Nuclide: ⁶⁸ Ge/ ⁶⁸ Ga

Check for transportation

- package meets all requirements for dangerous goods according to ADR/RID and IATA for class 7
- package is packed and labeled by trained and certified packer according to ADR/RID and IATA for class 7
- release of shipment from addressee submitted (EZR internal only)

Other: we would like to inform you, that in case the generator will not be returned as described Eckert & Ziegler will reserve the right to charge for any additional costs related to necessary extra safety arrangements!

Date:		Name (block letters)	
		Signature	